

ELECTRIC HOSPITAL BED HOME CARE; KL-HOME 03



Function

- Back up & down

- Knee up & down
- Turn left
- Turn right
- Integral up & down
- Defecation
- Washing head

Available Functions

Backrest lifting angle 0-85°

Knee-rest lifting angle 0-45°

Turn left and right 0-30°

Hi-lo adjustment
between 480mm-750mm

Packing list

Bed Board :Strip cold-rolled steel*1 pc

Head & Foot Board :ABS plastic detachable *1 pair

Side Rails:Foldable aluminum alloy*1 pair

Head Washing Basin:ABS plastic*1 pc

Bed Pan:ABS plastic*1 pc

Independent Brake Mute Catersl*4 pcs

Infusion Pole:Stainless steel*1 pc

Dining Table:High quality ABS plastic*1 pc

Motor:Branded powered motor,3 year quality guarantee;

Technical

Overall Dimensions L2150mm*W980mm*H 480mm-700mm
Hospital Bed Frame Perforated power coated steel
Hospital Bed Weight 120kg
Hospital Bed Platform MS perforated platform with ventilation hole
Loading Capacity 240kg (500Ibs)
Electric Bed Voltage: 10V/60Hz or 220V/50Hz
Options **Medical mattress; Bedside table; Bedside cabinet**
Packaging Box L2000mm*W1000mm*H500mm



ZÈÒÁ~ ÁÓÙÒÁ UWÙÒÛÄÊ€ÁË^} ^ Á~ Á& [| } ^|ÄÜŠËÛË ÕWYÄÄ
JHG €ÀVÖË ÙÄÜË ÒÒ
V^|KÉHFI JJÌ Ì Ì GÄ HJÍ HI FJH
ÖæKÉHFÍ JHG Í ÄÄ [àKÉHÍ Ì Ì GFÌ F
•[&ã ^ãO •[&ã ^ã& {