



LS-B35/50/75/100L-II (Hand Round Fully Automatically)

Characteristics:

1. Fully stainless steel structure.
2. Hand wheel type of quick-open door structure.
3. Door safety lock system.
4. Automatically controlled by computer.
5. Digital display of working status, touch of key.
6. Over temperature & over pressure auto-protection.
7. Steam-water inner circulation system: No steam discharge, and the environment for sterilizing will be clean and dry.
8. Self-inflating type seal.
9. Safe protection of water lacking.
10. Possible to install drying system according request.
11. Automatically shut off with beep reminding after sterilization.
12. With two stainless steel sterilizing baskets.

Specification:

Model	LS-B35L-II	LS-B50L-II	LS-B75L-II	LS-B100L-II
Valid volume	35L(φ318*450mm)	50L(φ340*550mm)	75L(φ400*600mm)	100L(φ440*650mm)
working pressure	0.22MPa			0.14MPa
Working temperature	134°C			126°C
Max working pressure	0.23 Mpa			0.165 Mpa
Heat average	± 1°C			
Timer	0~99min or 0~99hour			
Adjustment of temperature	0~134°C			0 ~ 126°C
Power	2.5KW/AC220V 50Hz	3KW /AC220V 50Hz	4.5KW /AC220V 50Hz	
Overall dimension	600*410*1030(mm)	650*410*1140(mm)	690*470*1140(mm)	730*510*1270(mm)
Transport dimension	730*500*1170(mm)	730*500*1290(mm)	760*550*1290(mm)	820*600*1380(mm)
G/N. W:	96Kg/77Kg	102Kg/ 83Kg	120Kg/98Kg	135/110Kg



Contacts

Z.A.C du BOIS MOUSSAY , 10 Avenue du colonel ROL-TANGUY - 93240
STAINS FRANCE

Tel: +33149988672, +33953341938

Fax: +33149988672 - Mob: +33678621481

socimed@socimed.com