

# BRITH 5080

## Electronic Delivery Table

BRITH 5080 is designed to provide at-most comfort and ease during delivery and inspection processes with all integrated functions and accessories ranging from pedal controlled height adjustment, backrest, Trendelenburg and reverse Trendelenburg positions to adjustable knee support pads, support handles. Combining all attributes with the exception of mobility, makes BRITH 5080 a revolutionary element of a maternity ward.



- 3 motors
- Electronic adjustments of backrest, height, Trendelenburg and Reverse Trendenbug functions by foot control pedal
- Anti-microbial, water proof, flame retardant, air permeable mattress
- Removable footrest mattress
- Chrome plated handles
- Stainless steel waste bin
- Lockable, adjustable and removable side rails
- Removable and Adjustable Polyurethane Leg Supports
- 360° Swivel Castors with central braking mechanism and steer control
- Protective plastic bumpers
- Detachable head board
- Height adjustable IV Pole
- ABS covered sub-frame
- Stainless steel accessory rails
- Easily cleanable, compact laminate mattress platform
- Electrostatic powder coated metal frame
- IV pole sockets

### Optional Accessories:

- Accessory trolley
- Battery Back-up
- Hand control unit
- Lifting pole
- Manual CPR
- File holder

TECHNICAL DETAILS	
<b>External Dimensions:</b>	min: 80x155cm
	max: 80x187cm
<b>Bed Surface:</b>	60x168cm
<b>Height (without mattress):</b>	62-91cm
<b>Backrest Angle:</b>	0°-70°
<b>Trendelenburg Angle:</b>	0°-14°
<b>Reverse Trendelenburg Angle:</b>	0°-14°
<b>Castor Diameter:</b>	150mm
<b>Safe Working Load Capacity:</b>	240kg
<b>Weight:</b>	93kg
<b>Mattress Thickness(backrest and pelvic sections)</b>	8cm
<b>Footrest Mattress Thickness:</b>	12cm
<b>Underbed Clearance:</b>	17cm
ELECTRICAL DATA	
<b>Power In:</b>	max. 3,15A 100V – 240V ~ 50/60 Hz
<b>Electrical Shock Protection:</b>	Class I, Type B
<b>Liquid Ingress Protection:</b>	IPx66