

# PX-300HF

Mobile X-ray unit



<b>Power</b>	<b>30kW (B -Type)</b>	<b>32kW (C -Type)</b>
<b>max. mA/kV</b>	<b>300mA/125kV</b>	<b>320mA/125kV</b>
<b>Tube</b>	<b>E7239X (1.0/2.0mm, 140kHU)</b>	<b>E7239X (1.0/2.0mm,140kHU)</b>
<b>PSU</b>	<b>Battery Type</b>	<b>Capacitor Type</b>