

ULTRA 200 & 300



ULTRA 300



ULTRA 200

- ✓ **Free stop Tube Arm**
- ✓ **180° lateral moving Tube Arm**
- ✓ **360° rotational wheels**
with multi-function pedals
- ✓ **Automatic SID display**
- ✓ **Automatic motor collimator**
- ✓ **Overhead OP console with touch screen**
- ✓ **100mA (ULTRA 200) output**
- ✓ **300mA (ULTRA 300) output**

Technical Specifications

Category	ULTRA 200	ULTRA 300
Output power	5 kW	20 kW
kV range	40 ~ 110 kV (1 kV step)	40 ~ 125 kV (1 kV step)
mA range	25 ~ 100 mA	10 ~ 300 mA
Focal spot	1.8 mm	Small : 0.6 mm / Large : 1.3 mm
Overhead OP console	Touch screen panel (4.3" LCD) Display : kV, mAs, tube head rotation, SID, tube temperature, wait time(sec.), field size, selectable APR(128)	
Motorized collimator	Auto SID by ultrasonic sensor Motorized auto collimation Max. field size : 47 x 47 cm @ 100 cm Collimator light : White LED Illumination : 130 LUX @ 100 cm	
UPS	Optional (200 exposures)	N/A
Line voltage	220 VAC (±10%)	110 VAC / 220 VAC (±10%)
Weight	166 kg	250 kg

