


HF-525Plus

General Radiographic Table X-ray System



Options

- Ceiling suspension tube stand
- 6 way table upgrade
- LCD OP Console upgrade
- PSU (Power Supply Unit)
- Automatic Exposure Control (AEC)
- DAP function



Power	32kW	40kW	50kW		65kW
Max mA/kV	400mA/125kV	500mA/125kV	630mA/150kV	800mA/150kV	800mA/150kV
Max rating	400mA @ 80kV	500mA @ 80kV	630mA @ 75kV	800mA @ 62kV	800mA@80kV
Tube	E7239X (1.0/2.0mm,140kHU)		E7252X (0.6/1.2mm,300kHU)		