



No.	Product Name	Parameters	Quantity (pc/pcs)	Image
1	Medical Monitor (BT240S)	24" HD Medical grade screen, Dimension ratio: 16:9 Input signal: DVI, VGA, CVBS Power supply:AC100~240V 50/60HZ Power:40VA Horizontal Resolution:1920*1080 Anti-jamming processing	1	
2	Medical Endoscope FHD Camera 1080P (SY-GW1000C-N)	Effective Pixels:1920*1080; Progressive scanning; Full hd camera 1080P; Anti-jamming processing; Horizontal Resolution:1080P(1200lines) Video Output: HDMI*4, SDI*1, BNC*1, VGA*1 Manual/ auto white balance ; SNR (Signal to Noise Ratio): ≥50dB Digital zoom: 0-60X ; Brightness: ODS adjusting Application: Laparoscopy, ENT, hysteroscopy, arthroscopy, urology, etc.	1	
3	HD Camera HHead	5 buttons, Freeze, White balance	1	
4	Optical Coupler	28mm Waterproof	1	
5	Led cold Light Source (SY-GW1000L-N)	LED light Power : 120watt Color temperature: ≥6500k Light image temperature: ≤80°C Digital display with cooling fan Working life is more than 20,000hours With timer function Application: Laparoscopy, ENT, hysteroscopy, arthroscopy, urology, etc.	1	
6	Light Guide Cable	2.5meter, standard diameter	1	
7	Co2 Insufflators (SY-GWQ300)	Rated power: 40VA; Range of flow:10-30L; Range of pressure: 5-30mmHg; Noise: ≤50dB; Power: 100-240V, 50/60Hz; Application: laparoscopy surgery	1	

8	Electrosurgical unit (SY-DP350)	<p>Adopting Automatic Power (P.A.C)system to maintain stable power when cutting tissues with different impedance.</p> <p>Three output ports:one for monopolar cut,one for monopolar fulguration and one for bipolar coagulation.</p> <p>Six working mode:regular cut .blended cut I,blended cut II,single spot coagulation,spurt coagulation amd bipolar coagulation.</p> <p>Independent power ofr electrical cut,electrical coagulation and bipolar coagulation</p> <p>Suitable for endoscopic surgery.</p> <p>Maximal out-put power: less than 400W;</p> <p>Six output modes</p>	1	
9	Trolley (SY-BT716)	<p>670*610*1310mm;</p> <p>Material: metal;</p> <p>Suitable robust and corrosion resistant trolley with five-storey shelves with suitable distance between the shelves.</p> <p>The trolley is provided with castors and two of the castors are lockable.</p>	1	
10	Endoscopic pump (SY-GWC200)	<p>Washing and suction the effluent or fragment in the enterocoelia surgery;</p> <p>Rate power: 50VA;</p> <p>Range of pos-press: 10~50kPa;</p> <p>Range of neg-press: -50~-10kPa ;</p> <p>Power: 100-240V, 50/60Hz;</p> <p>Application: laparoscopy surgery</p>	1	