


Ureter Stents

- Tumor Stents

MRI/EMG compatible

Made of Polyurethane
Tecoflex® USA
Laying time up to 30 days



Tumor Stent

Set consists of stent, pusher, suture and clamp.
Stent both sides open

Tumor Stents

TSS-624	6	24	0.035"
TSS-626	6	26	0.035"
TSS-628	6	28	0.035"
TSS-724	7	24	0.038"
TSS-726	7	26	0.038"
TSS-728	7	28	0.038"

Set with stainless steel guidewire add +G (for example EP-612-24+G).
Set with nitinol guidewire add +GTHS (for example EP-612-24+GTHS).

Class IIb device – for single use only