

# NASO-PHARYNGO-LARYNGOSCOPE

## OP-30

HIGH RESOLUTION IN THE SMALLEST DIAMETER

- Up to 18.000 pixels
- Detailed vision
- High resistance




# NASO-PHARYNGO-LARYNGOSCOPE OP-30



FIBROLUX DUO  
LIGHT SOURCE



FIBROLUX LED  
LIGHT SOURCE



FIBER OPTIC CABLE



FIBERSCOPE HOLDER



LEAKAGE TEST DEVICE



TWO-HANDED GRIP

	OP-30
Image guide system	18.000 pixels
Outer diameter	3,4mm
Flexible part diameter	3,2mm
Depth of field	2,5-50mm
Field of view	85°
Deflection angle	130°up/down
Direction of view	0°
Working length	300mm

- High definition images (up to 18.000 pixels).
- Focusing adjustment.
- ETO valve.
- Soakable in disinfectant liquid.
- Light guide connector for Storz and ACMI (optional Wolf).
- Delivered in a protective carrying case.