

Improving Patients Management *Improving Patients Experience*

HUBI-QUAN PRO[®]

One system for multiple markers

Make critical decisions on time with HUBI-QUAN PRO

HUBI-QUAN PRO[®] is a system to carry multiple biomarkers including Cardiac, Hormone, Inflammation, Tumor that provides the reliable test results to make rapid, cost-effective decisions for better patients treatment at the **point-of-care (POC)**.



- ✓ Simple test procedure
- ✓ Fast result in 5~15 minutes
- ✓ Multiple assays on one analyzer
- ✓ Auto-Calibration & Specimen recognition
- ✓ LIS/HIS Connectivity

User Instructions




Start run test



Drop the specimen on cartridge



Analyzing the cartridge



Result in 5~15 min.

Technical Information

Specifications

Dimension	W 255 x D 178 x H 92.7 mm, 1.35Kg	Printer	Built-in thermal, 203 dpi, 2 inch width
Processor	AMD PCS 498MHZ	Communication Port	USB(2ports), Serial port(RS-232), PS2, LAN
Memory	4GB	Voltage/Frequency	100~240V, 50~60HZ
Display	4.5 TFT LCD	Power Adaptor	Medical Power Supply, DC 12V output

Ordering Information

Analyzer

Category	Cat. No.	Product	Test time	Warranty Period
POCT Analyzer	AHQ-8001	HUBI-QUAN Pro	5 ~ 15 minutes	1 year

Test Kits

Category	Product	Cat. No.	Specimen	Shelf life	Category	Product	Cat. No.	Specimen	Shelf life
Myocardial Marker	HUBI Troponin I	ACTI-8025	Whole Blood, Plasma	18 months	Infectious Disease	HUBI CRP	ACRP-8025	Whole Blood (Capillary or Venous)	18 months
	HUBI CK-MB	ACCK-8025	Whole Blood, Plasma	12 months		HUBI PCT	APCT-8025	Whole Blood, Serum, Plasma	14 months
	HUBI BNP	ABNP-8025	Whole Blood, Plasma	18 months		HUBI IL-6	AIL-6-8025	Whole Blood, Serum, Plasma	18 months
	HUBI D-Dimer	ADIM-8025	Whole Blood, Plasma	18 months	Hormone	HUBI hCG plus	ANPP-8025	Whole Blood	24 months
	HUBI FABP	AFABP-8025	Whole Blood, Plasma	18 months		HUBI LH	AOV-8025	Venous blood	12 months
	HUBI DUO (TnI/CK-MB)	ACDC-8025	Whole Blood, Plasma	18 months		HUBI FSH	AME-8025	Venous blood	12 months
	HUBI Cardiac 3 in 1	ACTM-8025	Whole Blood, Plasma	18 months		HUBI TSH	ATSH-8025	Whole Blood, Serum, Plasma	12 months
	HUBI 3 in 1(B)	ACTCB-8025	Whole Blood, Plasma	18 months					
Tumor Marker	HUBI BPHScreen	ABPH-8025	Whole Blood (Capillary or Venous)	24 months					
	HUBI TESTO (coming soon)	ATST-8025	Whole Blood, Serum, Plasma	12 months					

Storage: 2-8 °C (except CRP/hCG plus: 2-30 °C)