



ENT Treatment Tables

Specialty Treatment Tables



SN-805 - For small areas
3 spray guns



SN-804 - Two physician work station
8 spray guns (4 on each side)

Features	SN-801	SN-802	SN-803	SN-804	SN-805
Over All Dimensions (LxWxH)	100x50x86 cm	100x50x86 cm	130x60x83 cm	120x70x86 cm	76x46x86 cm
Housing Stainless Steel #304	Yes	Yes	Yes	Yes	Yes
Spray Gun	4 pcs	5 pcs	5 pcs	8 pcs	3 pcs
Air Compressor Type (Oil-less, Maint. free)	Pressure Tank	Pressure Tank	Pressure Tank	Pressure Tank	Pressure Collector
Suction Pump (Oil-less, Maint. free)	Yes	Yes	Yes	Yes	Yes
Moveable Free Arm Light	Yes	Yes	Yes	Yes	Yes
Hand Control Eustachian Tube Catheter	1 pc	1 pc	1 pc	2 pcs	None
Power Source	110V, or 220V, or 240V 50/60 Hz				
Power Consumption	110V / 10A, or 220V / 5A, or 240V / 5A		110V / 12A, or 220V / 6A, or 240V / 3A, or 240V / 3A		
Fog-free Device	Yes	Yes	Yes	Yes - 2 sets	No
Eustachian Tube Catheter Device	Optional	Optional	Yes	Yes	No
Mirror Fog-free Device With Timer - 350W	Yes	Yes	Yes	Yes	No
Accessories					
Medicine Bottle (standard)	16 pcs	16 pcs	16 pcs	16 pcs	8 pcs
Moveable Instrument Plate (standard)	1 pc	1 pc	1 pc	1 pc	1 pc
Moveable Instrument Cup (standard)	2 pcs	2 pcs	2 pcs	2 pcs	1 pcs
Cotton Case (standard)	1 pc	1 pc	3 pcs	2 pcs	2 pcs
Electric Cauterizer Set (optional)	Available	Available	Available	Available	N/A

ENT Treatment Chair

Manual
seat back adjustment



Electric Hydraulic Pump
seat back adjustment



Specifications

Model	SN-602	SN-601
Dimensions (LxWxH)	78 x 62 x 56-84 cm	
Back Seat Angle	20° ~ 90°	
Back Seat Adjustment	Hydraulic Pump	Manual
Mattress	PVC foam covered with PVC plastic	
Elevation	56~84 cm	
Elevation Adjustment	Hydraulic Pump	
Headrest Adjustment	Manual	
Power Consumption	110V or 220V or 240V 50/60Hz 110V / 7A, 220V / 4A, 240V / 4A	

Features

1. Treatment chair controlled manually (SN-601) or by electric hydraulic system (SN-602).
2. Seat back, head rest and seat elevation are all adjustable.