

Tube strippers

Manual Tube stripper

Safe and easy tube stripping



The TS-Series Blood Tubing Stripper offers a comprehensive range of performance enhanced capabilities.

It is specifically designed to strip, crimp and cut tubing in the preparation of blood segments from the blood collection and processing procedures. The blood in the tubing segments is stripped into collection bags and mechanically sealed with the aluminum clips. Using the side jaws of the Hand Stripper, the segments are further separated by cutting between the seals.



Product Features

- High quality aluminum alloy
- Anti-rust, light weight and sturdy
- Ease of cleaning by wiping or soaking in IPA
- Ergonomic design for good grip and stripping
- Ease fatigue over prolong use
- Self-centering roller system keeps tubing in position
- Smooth traction for effective stripping of blood components
- Internal spring system keeps handles in open position
- Minimize entanglement with blood tubing



TS - 003

TS - 002

Technical data	TS - 002	TS - 003
Base Material	Aluminum Alloy & Stainless Steel	
Rollers	Special Composite Steel, Non-Rust	
Cutters	Not Applicable	Stainless Steel
Features	Self-Centering Rollers	
Cutter	No	Yes
Net Weight (g) +/-5%	68	70