

MEDICAL OXYGEN CONCENTRATOR

MODEL: JAY-10



DUAL FLOWS
USA imported humidifier
6PSI pressure release



Innovative "cooling"
technology to improve
reliability and life of the unit



Big LCD display
24H*365D work continuously



Optional

1. Nebulizer (Atomization)
2. Low purity alarm/Low flow alarm
3. High temperature alarm
4. Pulse oximeter
5. Maintenance alarm

Specification JAY-10

Flow Rate	0-10LPM
Purity	93% (±3%)
Outlet Pressure	0.04-0.07MPA
Sound Level	≤50db
Power Consumption	≤500w
LCD Display	Switch times, operating pressure, present working time, accumulated working time, presetting time from 10mins to 40hours
Alarm	Power Failure alarm, High&low pressure alarm
Net Weight	23Kgs / 25 Kgs / 28 Kgs
Size	365×375×600mm