

Gynecology and Obstetrics Operating Table

YOUR IDEAL OPERATING TABLE SUPPLIER

WELL SOLD FOR 35 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL SURGERY

ET400A Gynecological Examination Table



Introduction

Endearing appearance, various colors of mattress available, optional pedal switch electrically and easily operate.

It has overcome the slowness matter of similar products.

CONTACT

Z.A.C du BOIS MOUSSAY, 10 Avenue du colonel ROL-TANGUY - 93240 STAINS FRANCE
Tel: +33149988672, +33953341938 Fax: +33149988672 - Mob: +33678621481
socimed@socimed.com

Gynecology and Obstetrics Operating Table

YOUR IDEAL OPERATING TABLE SUPPLIER

WELL SOLD FOR 35 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL SURGERY

ET400C Gynecological Examination Table



Main Technical Data

Parameters may have reasonable error.

ET400A	ET400C
Length and width for tabletop: $\geq 1800/\geq 540\text{mm}$	Length and width for tabletop: $\geq 1780/\geq 600\text{mm}$
Elevation(down/range): $\leq 680/\geq 270\text{mm}$	Elevation(down/range): $\leq 680/\geq 300\text{mm}$
Trendlenburg/anti-trendelenburg: $\geq 17^\circ/\geq 14^\circ$	Trendlenburg/anti-trendelenburg: $\geq 20^\circ/\geq 8^\circ$
Back plate(up/down): $\geq 60^\circ/\geq 10^\circ$	Back plate(up/down): $\geq 60^\circ/\geq 7^\circ$
Power: 220V, 50Hz, 1.0 kW	Power: 220V, 50Hz, 1.0 kW

CONTACT

Z.A.C du BOIS MOUSSAY, 10 Avenue du colonel ROL-TANGUY - 93240 STAINS FRANCE
Tel: +33149988672, +33953341938 Fax: +33149988672 - Mob: +33678621481
socimed@socimed.com