



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

ET100 Electro-hydraulic Operating Table

Developments made in the conjunction with doctors: this is one of the basic principles which has made TECHARTMED the world market leader for operating tables. Close cooperation and years of experience have resulted in tried and tested products which facilitate the work in the operating room and help to increase the safety and comfort of the patient.



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery



CONTACT
Z.A.C du BOIS MOUSSAY , 10 Avenue du
colonel ROL-TANGUY - 93240 STAINS FRANCE
Tel: +33149988672, +33953341938
Fax: +33149988672 - Mob: +33678621481
socimed@socimed.com



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

ET100

Electro-hydraulic Operating Table



Your Ideal Operating Table Supplier

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

ET100

Electro-hydraulic Operating Table



Details

- 1, Leg plate: Imported gas spring for easy controls.
- 2, Two joint head plate.
- 3, Built-in kidney bridge.
- 4, Extreme minimum height 520mm.



TABLE HIGHEST HEIGHT



TABLE LOWEST HEIGHT



Product Description

ET100 as an advanced automatic operating table, it can meet various requirements of operations. Tabletop is X-ray available. Extremely low tabletop, especially suitable for cerebral surgery. Engineered for increasingly complex surgical scenes, ET100 comes with intelligent automatic control system satisfying various demands of general surgeries. Detachable head plate ensure the table can be use with the head stock.

- 520mm minimum and 820mm maximum heights for comfortable working, both seating and standing.
- Memory pad with antistatic, waterproof design.
- Optional battery power for 50~80 operations.

Well Sold For 40 Years Worldwide
Meet The Basic Demands Of General Surgery

Your Ideal Operating Table Supplier

| Model | Site | Length | Width | Tabletop lowest height | Elevation range |
|-------|------|---------|--------|------------------------|-----------------|
| ET100 | | ≥2085mm | ≥550mm | ≤520mm | 300mm |

PARAMETER ET100 electro-hydraulic operating table

| Movements | | Model: ET100 |
|-------------------------------|-----------------------|--------------|
| Tabletop trendelenburg | Reverse Trendelenburg | ≥25° |
| | Trendelenburg | ≥25° |
| Lateral Tilt (Left And Right) | | ≥15° |
| Head plate | Up | ≥40° |
| | Down | ≥90° |
| Back plate | Up | ≥80° |
| | Down | ≥20° |
| Leg plate | Up | ≥15° |
| | Down | ≥90° |
| Outward | | ≥90° |
| Kidney Bridge Elevation | | 110 mm |

STANDARD ACCESSORIES

| | | | | | | | |
|----------------------|---------|-------------------|--------|--------------|----------|-------------|---------|
| Anaesthesia screen | 1piece | Body support | 1 pair | Leg support | 1 pair | Hand remote | 1 piece |
| Kidney bridge handle | 1 piece | Arm support | 1 pair | Foot support | 1 pair | Mattress | 1 set |
| Shoulder support | 1 pair | Long fixing clamp | 1 pair | Fixing clamp | 8 pieces | Power Line | 1piece |

