

ET300C

Electric Operating Table

ET300C can meet various requirements from all departments surgery. Extra wide tabletop, long horizontal sliding that can be suitable for both X-ray and C-arm use. Adopted micro touch remote control that enables flexible and smooth adjustments on back plate and seat plate. With built-in kidney bridge. High automation, low noise, high reliability. Key parts adopted imported ones, can be regarded as an ideal electric operating table.



STANDARD ACCESSORIES

Anesthetic Screen	1 Piece
Kidney Bridge Handle	1 Piece
Shoulder Support	1 Pair
Body Support	1 Pair
Arm Support	1 Pair
Long Fixing Clamp	1 Pair
Leg Support	1 Pair
Foot Plate	1 Pair
Fixing Clamp	8 Pieces
Hand Remote	1 Piece
Mattress	1 Set
Power Line	1 Piece



CONTACT

Z.A.C du BOIS MOUSSAY , 10 Avenue du colonel ROL-
TANGUY - 93240 STAINS FRANCE
Tel: +33149988672, +33953341938
Fax: +33149988672 - Mob: +33678621481
socimed@socimed.com



YOUR IDEAL OPERATING TABLE SUPPLIER
MEET THE VARIOUS REQUIREMENTS OF GENERAL SURGERY

ET300C
Electric Operating Table





ET300C

Electric Operating Table

Optional 1006
Orthopedic Traction



Optional 1006
Tibia Bone-Nailing Position



Optional Lateral Arm Support



Optional Waste
Basin for Gynecology



Anorectum Operation Position



E.N.T. Operation Position



Built-in Kindey Bridge



Tabletop Sliding



Optional Head Rest



Elevation Up/Down



Lateral Tilt Left/Right



Trendelenburg/Reverse
Trendelenburg



Item	Model	ET300C
Tabletop Length		2070mm
Tabletop Width		550mm
Tabletop Minimum Height		700mm
Tabletop Maximum Height		1000mm
Trendelenburg/Reverse Trendelenburg		25°/25°
Lateral Tilt		15°/15°
Head Rest(Up/Down)		45°/90°
Foot Plate(Up/Down)		15°/90°
Foot Plate Horizontal Positioning		90°
Back Plate(up/down)		75°/20°
Kidney Bridge Elevation		110mm
Horizontal Sliding		300mm
Voltage		220V/110V
Frenquency		50Hz/60Hz
Power Capacity		1.0Kw
Battery		Yes
Mattress		Memory Mattress
Accessory		PU By Moulding
Main Material		304 Stainless Steel
Maximum Load Capacity		200KG
Warranty		2 Years

Movements Diagram of ET300C

