

## ET400B

Electric Table for Gynaecology and Obstetrics

ET400B is the most latest electric table that can be used for Obstetrics birthing , Gynecological operation and examination. The tabletop elevation, trendelenburg/reverse trendelenburg, back plate all controlled by both micro touch hand remote control and foot switch. Imported electric motors enable flexible, smooth movements, low noise, high reliability.

All base and column cover is made of 304 medical grade stainless steel. Seamless mattress can be colorfully customized, easy to clean and sterilize.



MEET THE VARIOUS  
REQUIREMENTS OF GYNAECOLOGY AND OBSTETRICS



YOUR IDEAL OPERATING TABLE SUPPLIER  
MEET THE VARIOUS REQUIREMENTS OF GYNAECOLOGY AND OBSTETRICS

## ET400B

Electric Table for Gynaecology and Obstetrics



### CONTACT

Z.A.C du BOIS MOUSSAY , 10 Avenue du  
colonel ROL-TANGUY - 93240 STAINS FRANCE  
Tel: +33149988672, +33953341938  
Fax: +33149988672 - Mob: +33678621481  
socimed@socimed.com



**ET400B**  
Electric Table for Gynaecology and Obstetrics

● Examination Position



● Birthing Position




**OPTIONAL ACCESSORY**

Examination Lamp

Various Colors For Mattress



**ET400B**  
Electric Table for Gynaecology and Obstetrics  
TABLETOP MOVEMENTS



**STANDARD ACCESSORIES**

|              |         |
|--------------|---------|
| Mattress     | 1 Set   |
| Arm Support  | 1 Pair  |
| Handle       | 1 Pair  |
| Fixing Clamp | 1 Pair  |
| Knee Crutch  | 1 Pair  |
| Pedal        | 1 Pair  |
| Leg Support  | 1 Piece |
| Waste Basin  | 1 Piece |
| Hand Remote  | 1 Piece |
| Foot Switch  | 1 Piece |

**MAIN SPECIFICATION**

| <b>ET400B</b> Electric Table for Gynaecology and Obstetrics |                 |
|-------------------------------------------------------------|-----------------|
| Tabletop Length/Width                                       | 1780/600mm      |
| Tabletop Up/Down                                            | 880/680mm       |
| Tabletop Trendelenburg/Reverse Trendelenburg                | 8°/20°          |
| Back Plate Up/Down                                          | 45°/7°          |
| Voltage/Frequency/Capacity                                  | 220V/50Hz/1.0kw |
| Emergency Battery                                           | Optional        |

**MOVEMENTS DIAGRAM OF ET400B**

