

## MT100 Manual Hydraulic Operating Table

### YOUR IDEAL OPERATING TABLE SUPPLIER

WELL SOLD FOR 35 YEARS WORLDWIDE  
MEET THE BASIC DEMANDS OF GENERAL SURGERY

MT100 Manual Hydraulic Operating Table



### Introduction

All mechanically operate, hydraulic lift. Back and leg plate are controlled by gas spring. Worm gear drive adjust the lateral tilt (left/right), trendelenburg and anti-trendelenburg movements. Simple structure and easily operate. The "T-shape" base design meets the ergonomics principle, makes enough free space for the doctors and very suitable for general surgery.

### Main Technical Data

Parameters may have reasonable error...

Item	MT100
Length and width for tabletop	≥2100/≥490
Elevation(down/range)	≤750/≥300mm
Trendlenburg/anti-trendelenburg	≥20°/≥30°
Lateral tilt(left/right)	≥20°/≥20°
Head plate(up/down)	≥45°/≥90°
Leg plate(down)	≥90°
Back plate(up/down)	≥75°/≥20°
Kidney bridge up	≥115mm

## MT300 Series Manual Hydraulic Operating Table

### YOUR IDEAL OPERATING TABLE SUPPLIER

WELL SOLD FOR 35 YEARS WORLDWIDE  
MEET THE BASIC DEMANDS OF GENERAL SURGERY

MT300A Manual  
Hydraulic Operating Table



MT300D Manual  
Hydraulic Operating Table



### Introduction

MT300A is widely used in chest, Abdominal surgery, ENT, gynaecology and obstetrics, urology and orthopedics etc.

Hydraulic lift by foot pedal, head operated movements.

The base and column cover all made by premium 304 stainless steel.

### Main Technical Data

Parameters may have reasonable error...

Item	MT300A
Length and width for tabletop	≥2080/≥480mm
Elevation(down/range)	≤790/≥240mm
Trendlenburg/anti-trendelenburg	≥25°/≥25°
Lateral tilt(left/right)	≥18°/≥18°
Head plate(up/down)	≥45°/≥90°
Leg plate(down/outward)	≥90°/0°
Back plate(up/down)	≥90°/≥10°
Kidney bridge up	≥120mm