



digital

HY5510

Laptop Ultrasonic Diagnostic System



socimed

VENTE DE CONSOMMABLES ET D'APPAREILS MÉDICAUX
Tel. +331 49 98 86 72 - socimed@socimed.com





HY5510

Laptop Ultrasonic Diagnostic System



Feature:


- 10 inch high definition color monitor
- 8 stage TGC, gain adjustment
- AC and DC power supply, adaptor
- Large capacity lithium battery
- Clinical measurement software
- Wide- band multi frequency probes
- Automatic probe recognition

Standard Configuration:

- Main Unit, 3.5Mhz Convex probe

Options:

- USB ports, printer function, 6.5Mhz transvaginal probe, 7.5MHz Linear probe, 7.5Mhz transrectal probe, etc.



socimed

VENTE DE CONSOMMABLES ET D'APPAREILS MÉDICAUX
Tel. +331 49 98 86 72 - socimed@socimed.com

