


FICHE PRODUIT



Plaque DCP auto compressive pour vis 3,5 mm 2 à 12 trous

Référence : SOC-040000

Description du produit

Plaque DCP auto compressive pour vis 3,5 mm ; 2 à 12 trous

- Nombre de trous 2 à 12
- Largeur plaque 10 mm
- Epaisseur plaque 3,5 mm
- Ecartement entre les trous 12 mm et 16 mm pour l'écartement médian

Conformité et garantie



Garantie

Garantie 2 ans

Informations réglementaires : Dispositif médical : vérifier la notice, l'étiquetage et la déclaration UE de conformité avant utilisation.

SOCIMED SARL

9, Chaussée Jules César (Bâtiment 01 Porte 106)
95520 OSNY FRANCE
Tél. : +33(0)1 55 93 12 57
fr.quotations@socimed.com

Page produit

<https://www.socimed.com/plaques-de-verrouillage-orthopedique/13806-plaque-dcp-auto-compressive-pour-vis-3-5-mm-2-a-12-trous.html>



Déclinaisons

Désignation	Référence :	EAN13 :
trous_plaque_coudee - 2 trous	SOC-0400021	
trous_plaque_coudee - 3 trous	SOC-040003	
trous_plaque_coudee - 4 trous	SOC-040004	
trous_plaque_coudee - 5 trous	SOC-040005	
trous_plaque_coudee - 6 trous	SOC-040006	
trous_plaque_coudee - 7 trous	SOC-040007	
trous_plaque_coudee - 8 trous	SOC-040008	
trous_plaque_coudee - 9 trous	SOC-040009	
trous_plaque_coudee - 10 trous	SOC-040010	
trous_plaque_coudee - 12 trous	SOC-040012	

SOCIMED SARL

9, Chaussée Jules César (Bâtiment 01 Porte 106)
 95520 OSNY FRANCE
 Tél. : +33(0)1 55 93 12 57
fr.quotations@socimed.com